



Per una riproduzione esatta dei colori è consigliabile fare riferimento ad una cartella colori RAL originale.
For a precise reproduction, we suggest you to make reference to an original copy of RAL color chart.
Para una reproducción exacta de colores, es aconsejable consultar una tabla de colores RAL original.

Colorazioni disponibili
Colors available
Colores disponibles



Cradle adjustable height on castors with reverse-Trend

Height adjustable cot through gas spring control. "X" shaped basement in order to increase the stability and to place the cot as closer as possible to the hospital beds.

Fixed to the central column equipped with twins wheels, two with brake.

Column made of extruded aluminium.

Swinging support system with absolute precision adjustment that allows reverse-Trendelenburg of 13°.

Basin made of non-toxic material, transparent and thermoshaped, provided with place for label and holes at the bottom for air circulation.

Can be provided with box for clothes or personal belongings.

Accessory:



SZ0098-COVER DOME MADE OF PETG FOR SL2110/SL2100 CRIBS FOR INFANTS



SZ0095-LOWER BASKET FOR SL2110 CRIB FOR INFANT