



Double side cabinet with 1 upper drawer, open space and one lower big drawer

COPERNICO HC8120 and HC8125 are equipped with one or two upper drawers, total extraction and 3kg capacity. The capacious lower drawer with total extraction it is easy to access by the patient with 7kg capacity and system against overturning. The upper shelf is shaped in order to contain liquids, easy to clean and disinfect while the bottom is provided with 4 holes to grant the drainage of liquids and the airing. Back equipped with towel support and bottle support (n° 2 bottles of 1,5lt each). In the middle there is the lock for the first drawer of the cabinet which become a safe place for personal belongings.

The handles are shaped in such way that can be used from both sides.

Accessory:



HZ0020-LOCK FOR FIRST DRAWER OF COPERNICO BEDSIDE CABINETS